

Product Datasheet

Spike glycoprotein antibody orb863907

Description	Spike glycoprotein antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human coronavirus 229E Spike glycoprotein(S),partial (785-873aa)
Target	S
Form/Appearance	Liquid
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P15423
Expiration Date	12 months from date of receipt.