

Product Datasheet

Spike glycoprotein antibody orb864028

Description	Rabbit polyclonal antibody to Spike glycoprotein. All of our unique antibodies are supplied with the following: 200 µg recombinant immunogen protein/peptide (Positive Control). 10mg Rabbit polyclonal antibody purified by Protein A/G.
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Bat coronavirus 133/2005 S protein
Target	S
Form/Appearance	Liquid
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q0Q4F2
Expiration Date	12 months from date of receipt.