

## Product Datasheet

### Nucleoprotein antibody orb864184

|                            |                                                                                                                                                                                                                                      |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to Nucleoprotein. All of our unique antibodies are supplied with the following: 200 µg recombinant immunogen protein/peptide (Positive Control). 10mg Rabbit polyclonal antibody purified by Protein A/G. |
| <b>Species/Host</b>        | Rabbit                                                                                                                                                                                                                               |
| <b>Reactivity</b>          | Virus                                                                                                                                                                                                                                |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                                                                         |
| <b>Tested Applications</b> | ELISA, WB                                                                                                                                                                                                                            |
| <b>Immunogen</b>           | Recombinant Bovine coronavirus N protein                                                                                                                                                                                             |
| <b>Target</b>              | N                                                                                                                                                                                                                                    |
| <b>Form/Appearance</b>     | Liquid                                                                                                                                                                                                                               |
| <b>Preservatives</b>       | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol                                                                                                                                                                          |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                                                                                                                                                        |
| <b>Note</b>                | For research use only.                                                                                                                                                                                                               |
| <b>Isotype</b>             | IgG                                                                                                                                                                                                                                  |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                                                                           |
| <b>Uniprot ID</b>          | <b>Q9QAR8</b>                                                                                                                                                                                                                        |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                                                                                      |