

Product Datasheet

Novel Coronavirus Nucleoprotein (N) antibody orb783580

Description	Mouse monoclonal antibody to Novel Coronavirus Nucleoprotein (N)
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, GICA
Immunogen	Recombinant Human Novel Coronavirus Nucleoprotein (N) (1-419aa) (CSB-EP3325GMY)
Target	N
Form/Appearance	Liquid
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only.
Isotype	Mouse scFv fusion with human IgG1 Fc
Clonality	Monoclonal
Clone ID	4D2
Uniprot ID	P0DTC9
Dilution Range	ELISA:1:10000-1:50000, LF:1:500-1:5000
Expiration Date	12 months from date of receipt.