

Product Datasheet

MARK2 (phospho-Thr596) antibody (Biotin) orb782919

Description	MARK2 (phospho-Thr596) antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Predicted Reactivity	Bovine, Canine, Equine, Guinea pig, Porcine, Sheep
Conjugation	Biotin
Tested Applications	ELISA, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human MARK2 around the phosphorylation site of Ser595
Target	phospho-MARK2 (Ser595)
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q7KZI7
Dilution Range	WB=1:500-2000,IHC-P=1:100-500,IHC-F=1:100-500,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.