

Product Datasheet

RBCK1 antibody orb783098

Description	RBCK1 antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant human RBCK1 (UBCE7IP3).
Target	RBCK1
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG1,k
Clonality	Monoclonal
MW	56 kDa
Clone ID	R5E2
Uniprot ID	Q9BYM8
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.