


Product Datasheet

KIF5B antibody orb783194

Description	KIF5B antibody
Species/Host	Rabbit
Reactivity	Rat
Predicted Reactivity	Bovine, Canine, Equine, Human, Mouse, Porcine, Sheep
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human KIF5B
Target	KIF5B
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	106 kDa
Uniprot ID	P33176
Dilution Range	IHC-P=1:400-800, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



IHC-P analysis of rat brain using KIF5B Polyclonal Antibody at 1:200 dilution.