

Product Datasheet

Lamin A antibody orb783219

Description	Lamin A antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant human Lamin A/C.
Target	Lamin A/C
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG1,k
Clonality	Monoclonal
MW	73 kDa
Clone ID	C3S4
Uniprot ID	P02545
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.