

Product Datasheet

ALDH6A1 antibody orb783407

Description	ALDH6A1 antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P, WB
Immunogen	Recombinant human ALDH6A1.
Target	ALDH6A1
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG1,IgK
Clonality	Monoclonal
MW	59 kDa
Clone ID	AL06
Uniprot ID	Q02252
Dilution Range	WB=1:500-1000, IHC-P=1:25, ICC=1:25, IF=1:25
Expiration Date	12 months from date of receipt.