

Product Datasheet

ALDH2 antibody orb783408

Description	ALDH2 antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-P, WB
Immunogen	Recombinant human ALDH2.
Target	ALDH2
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG1,Ig?
Clonality	Monoclonal
MW	54 kDa
Clone ID	A02C1
Uniprot ID	P05091
Dilution Range	WB=1:500, IHC-P=1:400-800, ICC=1:50, Flow-Cyt=1:50, IHC-F=1:400-800, IF=1:100-500
Expiration Date	12 months from date of receipt.