

Product Datasheet

C-RAF antibody (Biotin) orb783481

Description	C-RAF antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Predicted Reactivity	Canine, Equine, Gallus, Mouse, Porcine, Rat
Conjugation	Biotin
Tested Applications	ICC, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human RAF1
Target	RAF1
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P04049
Dilution Range	IHC-P=1:100-500,IHC-F=1:100-500,ICC=1:100-500
Expiration Date	12 months from date of receipt.