

Product Datasheet

Nrf2 (phospho-Ser40) antibody (Biotin) orb783483

Description	Nrf2 (phospho-Ser40) antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human, Mouse
Predicted Reactivity	Bovine, Canine, Gallus, Rat
Conjugation	Biotin
Tested Applications	ELISA, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Nrf2 around the phosphorylation site of Ser40
Target	phospho-Nrf2 (Ser40)
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q16236
Dilution Range	WB=1:500-2000,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.