

Product Datasheet

Growth Hormone Receptor antibody (Biotin) orb783497

Description	Growth Hormone Receptor antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Predicted Reactivity	Bovine, Canine, Gallus, Porcine, Rabbit, Sheep
Conjugation	Biotin
Tested Applications	ELISA, WB
Immunogen	KLH conjugated synthetic peptide derived from human GHR
Target	Growth hormone receptor
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P10912
Dilution Range	WB=1:500-2000,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.