

Product Datasheet

HSP90 Beta antibody (FITC) orb783499

Description	HSP90 Beta antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Predicted Reactivity	Bovine, Equine, Mouse, Porcine, Rat
Conjugation	FITC
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from human HSP90 beta
Target	HSP90 beta
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P08238
Dilution Range	Flow-Cyt=1?g/Test
Expiration Date	12 months from date of receipt.