

Product Datasheet

Anti-Human secretory IgA antibody {HRP} orb100834-HRP

Description	Rabbit Polyclonal to Human IgA.
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ICC, IF
Immunogen	Human IgA
Target	IgA
Form/Appearance	Liquid: PBS, 2 mM Sodium Azide
Concentration	100 µg (in 411 µl): 0.24 mg/ml
Storage	Aliquot and store at -20 °C or -80 °C, avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified
Dilution Range	To be determined by the end user
Expiration Date	12 months from date of receipt.