

Product Datasheet

GNPTAB antibody {HRP} orb157234-HRP

Description	Rabbit polyclonal antibody against GNPTAB conjugated to HRP.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Predicted Reactivity	Porcine
Conjugation	HRP
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from the region aa11-90 of human GNPTAB
Target	GNPTAB {HRP}
Form/Appearance	Liquid: PBS and 2mM Sodium Azide
Concentration	100 µg (in 411 µl): 0.24 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified by Protein A
MW	105 kDa
Uniprot ID	Q3T906 (Human)
Entrez	79158
Dilution Range	To be determined by the end user
Expiration Date	12 months from date of receipt.