

Product Datasheet

SV40 VP1 antibody {HRP} orb313227-HRP

Description	Rabbit polyclonal antibody against SV40 VP1 conjugated to HRP
Species/Host	Rabbit
Conjugation	HRP
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from the region aa261-340 of human SV40 VP1
Target	SV40 VP1 {HRP}
Form/Appearance	Liquid: PBS, 2 mM Sodium Azide
Concentration	100 µg (in 411 µl): 0.24 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified by Protein A
MW	40 kDa
Dilution Range	To be determined by the end user
Expiration Date	12 months from date of receipt.