

Product Datasheet

phospho-GABBR2 (Ser884) Polyclonal Antibody orb1088899

Description	phospho-GABBR2 (Ser884) Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Animal, Bovine, Canine, Guinea pig, Human, Mouse, Porcine, Rabbit, Rat, Sheep
Tested Applications	WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human GABBR2 around the phosphorylation site of Ser884:IN(p-S)PE
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only.
Clonality	Polyclonal
MW	104 kD
Uniprot ID	O75899
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.