

Product Datasheet

phospho-ERK1 + 2 (Thr183/Tyr185) Polyclonal Antibody orb1088867

Description	phospho-ERK1 + 2 (Thr183/Tyr185) Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Animal, Bovine, Canine, Gallus, Guinea pig, Human, Mouse, Rabbit, Rat
Tested Applications	ICC, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from mouse ERK1 around the phosphorylation site of Thr183/Tyr185:FL(p-T)E(p-Y)V
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only.
Clonality	Polyclonal
MW	42 kD
Uniprot ID	P63085
Dilution Range	WB=1:500-2000,ICC=1:100-500
Expiration Date	12 months from date of receipt.