

Product Datasheet

Neutrophil Elastase antibody {CF488A} orb500993-CF488A

Description	Rabbit polyclonal antibody against Neutrophil Elastase to CF488A.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	CF488A
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from the region aa111-190 of human Neutrophil Elastase
Target	Neutrophil Elastase {CF488A}
Form/Appearance	Liquid: PBS and 2mM Sodium Azide
Concentration	100 µg in 500 µl (0.20 mg/ml)
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified by Protein A
MW	26 kDa
Uniprot ID	P08246 (Human); Q3UP87 (Mouse)
Entrez	1991
Dilution Range	To be determined by end user
Expiration Date	12 months from date of receipt.