

Product Datasheet

Egr1 antibody {CF405M} orb10586-CF405M

Description	Rabbit polyclonal antibody against Egr1 conjugated to CF405M.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	CF405M
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived between 411-490 amino acids of human Egr1
Target	Egr1 {CF405M}
Form/Appearance	Liquid: PBS, 2 mM Sodium Azide, 30% glycerol
Concentration	100 µg in 500 µl (0.20 mg/ml)
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified by Protein A
MW	60 kDa
Uniprot ID	P18146 (Human); P08046 (Mouse)
Entrez	1958
Dilution Range	To be determined by end user
Expiration Date	12 months from date of receipt.

