
Product Datasheet

20x Tris-Buffered Saline with Tween 20, pH:7.6 orb1109222

Description

ClearBand 20x Tris-Buffered Saline with Tween 20 (TBST) is a space-saving balanced pH stabilizing salt solution used for especially western blotting and ELISA procedures. ClearBand TBST enables washing without disturbing antibody-antigen binding interactions. ClearBand 20x TBST should be diluted to 1x working solution. 1x formulation contains 0.1% Tween 20.

Note

For research use only.

Expiration Date

12 months from date of receipt.