

Product Datasheet

HoxC13 antibody {CF405M} orb157569-CF405M

Description	Rabbit polyclonal antibody against HoxC13 conjugated to CF405M.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Predicted Reactivity	Bovine, Canine, Equine, Gallus, Sheep
Conjugation	CF405M
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from the region aa241-320 of human HoxC13
Target	HoxC13 {CF405M}
Form/Appearance	Liquid: PBS, 2 mM Sodium Azide, 30% glycerol
Concentration	100 µg in 500 µl (0.20 mg/ml)
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified by Protein A
MW	35 kDa
Uniprot ID	P31276 (Human); P50207 (Mouse)
Entrez	3229
Dilution Range	To be determined by end user
Expiration Date	12 months from date of receipt.