

Product Datasheet

BLM antibody {HRP} orb155843-HRP

Description	Rabbit polyclonal antibody to BLM.
Reactivity	Human, Mouse, Rat
Predicted Reactivity	Canine, Equine, Bovine, Porcine, Sheep
Conjugation	HRP
Tested Applications	ELISA, ICC, IF, IHC-P, WB
Concentration	100 µg in 500 µl (0.20 mg/ml)
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Dilution Range	To be determined by the end user
Expiration Date	12 months from date of receipt.