

## Product Datasheet

### Mouse IgM Isotype Control for In Vivo - Low Endotoxin ICH2252

<b>Description</b>	Mouse IgM Isotype Control for In Vivo - Low Endotoxin
<b>Species/Host</b>	Mouse
<b>Concentration</b>	? 2.0 mg/ml
<b>Preservatives</b>	This Low Endotoxin Functional Formulation antibody is formulated in 0.01 M phosphate buffered saline (PBS) pH 7.2 - 7.4.
<b>Storage</b>	This Low Endotoxin Functional Formulation antibody is stable when stored at 2-8°C. Do not freeze.
<b>Note</b>	For research use only.
<b>Isotype</b>	IgM
<b>Clonality</b>	Monoclonal
<b>Purity</b>	> 95% by SDS-PAGE and HPLC
<b>Expiration Date</b>	12 months from date of receipt.
<b>Application Notes</b>	This isotype control antibody should be used at the same concentration as the primary antibody.