

Product Datasheet

Bulk anti-Human CD5 antibody (UCHT-2) ICH1174

Description	Bulk anti-Human CD5 antibody (UCHT-2)
Species/Host	Mouse
Reactivity	Human
Tested Applications	FC, IF, IHC-Fr, FA
Target	CD5
Concentration	? 5.0 mg/ml
Preservatives	This monoclonal antibody is aseptically packaged and formulated in 0.01 M phosphate buffered saline (150 mM NaCl) PBS pH 7.2 - 7.4 with no carrier protein, potassium, calcium or preservatives added.
Storage	This antibody is stable for at least one week when stored sterile at 2-8°C. For long term storage aseptically aliquot in working volumes without diluting and store at -80°C. Avoid Repeated Freeze Thaw Cycles.
Note	For research use only.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	?95% monomer by analytical SEC, > 95% by SDS Page
Clone ID	UCHT-2
Uniprot ID	P06127
Expiration Date	12 months from date of receipt.
Application Notes	The suggested concentration for this CD5 antibody for staining cells in flow cytometry is ? 1.0 ?g per 10 ⁶ cells in a volume of 100 ?l. Titration of the reagent is recommended for optimal performance for each application.