

Product Datasheet

Bulk anti-Human CD45RO antibody (UCHL-1) ICH1175

Description	Bulk anti-Human CD45RO antibody (UCHL-1)
Species/Host	Mouse
Reactivity	Human
Tested Applications	FC, IF, IHC-Fr, IHC-P, IP, WB
Target	CD45RO
Concentration	? 5.0 mg/ml
Preservatives	This monoclonal antibody is aseptically packaged and formulated in 0.01 M phosphate buffered saline (150 mM NaCl) PBS pH 7.2 - 7.4 with no carrier protein, potassium, calcium or preservatives added.
Storage	This antibody is stable for at least one week when stored sterile at 2-8°C. For long term storage aseptically aliquot in working volumes without diluting and store at -80°C. Avoid Repeated Freeze Thaw Cycles.
Note	For research use only.
Isotype	Mouse IgG2a
Clonality	Monoclonal
Purity	?95% monomer by analytical SEC, > 95% by SDS Page
Clone ID	UCHL-1
Uniprot ID	P08575
Expiration Date	12 months from date of receipt.
Application Notes	The suggested concentration for this CD45 antibody for staining cells in flow cytometry is ? 1.0 ?g per 10 ⁶ cells in a volume of 100 ?l. Titration of the reagent is recommended for optimal performance for each application.