

Product Datasheet

Bulk Human IgG1 Isotype Control (IB1) ICH2254

Description	Bulk Human IgG1 Isotype Control (IB1)
Species/Host	Human
Reactivity	Human
Tested Applications	FC, FA
Target	IgG1
Concentration	? 5.0 mg/ml
Preservatives	This monoclonal antibody is aseptically packaged and formulated in 0.01 M phosphate buffered saline (150 mM NaCl) PBS pH 7.2 - 7.4 with no carrier protein, potassium, calcium or preservatives added.
Storage	This antibody is stable for at least one week when stored sterile at 2-8°C. For long term storage aseptically aliquot in working volumes without diluting and store at -80°C. Avoid Repeated Freeze Thaw Cycles.
Note	For research use only.
Isotype	Human IgG1
Clonality	Monoclonal
Purity	?95% monomer by analytical SEC, > 95% by SDS Page
Clone ID	IB1
Expiration Date	12 months from date of receipt.