

Product Datasheet

Mouse Anti-Mouse PD-1 antibody (RMP1-14) ICH1182

Description	Mouse Anti-Mouse PD-1 Antibody (RMP1-14)
Species/Host	Mouse
Reactivity	Mouse
Tested Applications	Blocking, FC, WB, FA
Immunogen	Mouse PD-1 transfected BHK cells
Target	PD-1
Concentration	1.0 - 5.0 mg/ml
Preservatives	0.01 M phosphate buffered saline (PBS) pH 7.2, 150 mM NaCl with no carrier protein, potassium or preservatives added.
Storage	Mouse anti-mouse PD-1 antibody (RMP1-14) is stable when stored at 2-8°C. For longer term storage (> 3 months) it is recommended to store this antibody at -70°C.
Note	For research use only.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	> 95% by SDS-PAGE and HPLC
Clone ID	RMP1-14
Uniprot ID	Q02242
Expiration Date	12 months from date of receipt.
Application Notes	Each investigator should determine their own optimal working dilution for specific applications.