

## Product Datasheet

### Tim-3 Biosimilar Antibody ICH5079

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Sabatolimab Biosimilar - Research Grade  |
| <b>Reactivity</b>          | Human  |
| <b>Tested Applications</b> | FA   |
| <b>Endotoxins</b>          | ? 1.0 EU/mg as determined by the LAL method  |
| <b>Target</b>              | Tim-3  |
| <b>Concentration</b>       | 5mg/ml   |
| <b>Preservatives</b>       | Sterile, preservative-free, solution in PBS. rAlbumin and Azide free.  |
| <b>Storage</b>             | This biosimilar is stable when stored at 2-8°C. For longer term storage (> 3 months) it is recommended to store this antibody at -20°C or -70°C. |
| <b>Note</b>                | For research use only.   |
| <b>Isotype</b>             | Human IgG4 S228P-Kappa   |
| <b>Clonality</b>           | Monoclonal   |
| <b>Purity</b>              | > 95% by SDS-PAGE  |
| <b>Uniprot ID</b>          | <b>Q8TDQ0</b>  |
| <b>Expiration Date</b>     | 12 months from date of receipt.  |
| <b>Application Notes</b>   | This biosimilar is for research use only (RUO): it is not for diagnostic or therapeutic procedures and cannot be purchased by patients.          |