

Product Datasheet

TNFSF13B Biosimilar Antibody ICH5080

Description	Belimumab Biosimilar - Research Grade
Reactivity	Human
Tested Applications	FA
Endotoxins	? 1.0 EU/mg as determined by the LAL method
Target	TNFSF13B
Concentration	5mg/ml
Preservatives	Sterile, preservative-free, solution in PBS. rAlbumin and Azide free.
Storage	This biosimilar is stable when stored at 2-8°C. For longer term storage (> 3 months) it is recommended to store this antibody at -20°C or -70°C.
Note	For research use only.
Isotype	Human IgG1 E239D ,M241L-Human lambda3 K83R
Clonality	Monoclonal
Purity	> 95% by SDS-PAGE
Uniprot ID	Q9Y275
Expiration Date	12 months from date of receipt.
Application Notes	This biosimilar is for research use only (RUO): it is not for diagnostic or therapeutic procedures and cannot be purchased by patients.