

Product Datasheet

Human CYP2E1 ELISA Kit orb404261

Description	Human Cytochrome P450 2E1(CYP2E1) ELISA kit
Reactivity	Human
Range	43.75 pg/mL-2800 pg/mL
Sensitivity	10.94 pg/mL
Target	CYP2E1
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P05181
Expiration Date	Please enquire.
Application Notes	Detection Wavelength: 450 nm