

Product Datasheet

Human CDKN2A ELISA Kit orb404243

Description	Human Cyclin-dependent kinase inhibitor 2A, isoforms 1/2/3(CDKN2A) ELISA kit
Reactivity	Human
Range	18.75 pg/mL-1200 pg/mL
Sensitivity	4.69 pg/mL
Target	CDKN2A
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P42771
Expiration Date	Please enquire.
Application Notes	Detection Wavelength: 450 nm