

Product Datasheet

Human LPAR3 ELISA Kit orb440593

Description	Human Lysophosphatidic Acid Receptor 3 (LPAR3) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Sensitivity	0.105ng/mL
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9UBY5
NCBI	23566
Entrez	23566
Expiration Date	Please enquire.
Application Notes	tissue homogenates, cell lysates and other biological fluids.