

Product Datasheet

SEPT10 antibody orb100049

Description	SEPT10 antibody
Species/Host	Rabbit
Reactivity	Mouse
Predicted Reactivity	Equine, Gallus, Human, Porcine, Rabbit, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human SEPT10 (121-220/454 aa)
Target	SEPT10
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	53 kDa
Uniprot ID	Q9P0V9
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.