

Product Datasheet

Cytokeratin 18 antibody {CF488A} orb5871-CF488A

Description	Rabbit polyclonal antibody to Cytokeratin 18
Reactivity	Bovine, Equine, Human, Mouse, Porcine, Rabbit, Rat
Predicted Reactivity	Canine, Equine, Bovine, Porcine
Conjugation	CF488A
Tested Applications	ICC, IF
Concentration	100 µg in 500 µl (0.20 mg/ml)
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Dilution Range	To be determined by the end user
Expiration Date	12 months from date of receipt.