

Product Datasheet

NFKB1 (phospho-Ser932) antibody (PE) orb504636

Description	NFKB1 (phospho-Ser932) antibody (PE)
Species/Host	Rabbit
Reactivity	Human, Rat
Predicted Reactivity	Bovine, Canine, Gallus, Mouse, Porcine
Conjugation	PE
Tested Applications	FC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human NF KappaB p105 around the phosphorylation site of Thr931
Target	Phospho-NFKB1 (Thr931)
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P19838
Dilution Range	Flow-Cyt=1?g/Test,IF=1:100-500
Expiration Date	12 months from date of receipt.