

Product Datasheet

Goat Anti-Rat IgG H&L antibody orb868164

Description	Goat Anti-Rat IgG H&L antibody
Species/Host	Goat
Reactivity	Rat
Conjugation	Unconjugated
Immunogen	Native Rat IgG
Form/Appearance	Lyophilized powder
Preservatives	0.01M PBS (pH7.4)
Storage	Store at -20°C stable for 2 years (lyophilized powder). Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.