

Product Datasheet

Rabbit Anti-Avidin antibody orb868162

Description	Rabbit Anti-Avidin antibody
Species/Host	Rabbit
Reactivity	All
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-Fr, IHC-P, WB
Immunogen	Native Avidin
Form/Appearance	Lyophilized powder
Preservatives	0.01M PBS (pH7.4)
Storage	Store at -20°C stable for 2 years (lyophilized powder). Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Dilution Range	WB=1:1000-5000,IHC-P=1:100-500,IHC-F=1:200-1000,ELISA=1:2000-10000
Expiration Date	12 months from date of receipt.