

## Product Datasheet

### EPOR antibody (PE) orb124613

|                      |                                                                                  |
|----------------------|----------------------------------------------------------------------------------|
| Description          | EPOR antibody (PE)                                                               |
| Species/Host         | Rabbit                                                                           |
| Reactivity           | Human, Rat                                                                       |
| Predicted Reactivity | Bovine, Canine, Equine, Mouse                                                    |
| Conjugation          | PE                                                                               |
| Tested Applications  | IF                                                                               |
| Immunogen            | KLH conjugated synthetic peptide derived from human EPOR                         |
| Target               | EPOR                                                                             |
| Form/Appearance      | Liquid                                                                           |
| Concentration        | 1mg/ml                                                                           |
| Preservatives        | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.            |
| Storage              | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. |
| Note                 | For research use only.                                                           |
| Isotype              | IgG                                                                              |
| Clonality            | Polyclonal                                                                       |
| Uniprot ID           | <a href="#">P19235</a>                                                           |
| Dilution Range       | IF=1:100-500                                                                     |
| Expiration Date      | 12 months from date of receipt.                                                  |