

Product Datasheet

IL1 alpha Receptor antibody (FITC) orb466067

Description	IL1 alpha Receptor antibody (FITC)
Species/Host	Rabbit
Reactivity	Human, Mouse
Predicted Reactivity	Rat
Conjugation	FITC
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from human IL-1R1
Target	IL-1R1
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P14778
Dilution Range	Flow-Cyt=1?g /test
Expiration Date	12 months from date of receipt.