

Product Datasheet

Anti-LAR/PE orb498778

Description	Anti-LAR/PE
Species/Host	Rabbit
Reactivity	Human
Predicted Reactivity	Bovine, Canine, Equine, Mouse, Rabbit, Rat, Sheep
Conjugation	PE
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human PTPRF
Target	LAR
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	210kDa
Uniprot ID	P10586
Dilution Range	Flow-Cyt=1?g/Test,IF=1:100-500,ICC=1:100-500
Expiration Date	12 months from date of receipt.